

Dalyell Group Week 4

Exercise

yu.sun@sydney.edu.au

1. Perform some manipulation similar to what is done in the recording. You may need to Google the relevant functions. The numpy [documentation](#) for n-dimension arrays can be helpful.

In a Python interactive shell:

- a. Read the koala image `koala.tiff`;
 - b. Find the 75 percentile (`p75`) of all signal intensities (SI);
 - c. Find the pixels whose SI is above `p75` and halve their SI;
 - d. Save the image using the file name `koala_processed.tiff`;
 - e. Get the difference image (the SI difference of the `koala_processed.tiff` compared with the `koala.tiff`) and save it as `koala_diff.tiff`;
 - f. Rotate `koala_diff.tiff` 90 degree clockwise and scale it to 50% of the original size. Save the result as `koala_diff_rot.tiff`.
2. Medical images are typically stored in a format call DICOM. A DICOM file has a header which contains the meta-information (i.e. scan type, scan date, etc.). To read a DICOM image, you need to use specialised modules, such as `pydicom`.

Install the module `pydicom`. Refer to this [guide](#) about how to install third-party Python modules. Typically you should get it installed in the terminal (cmd in Windows) using

```
pip install pydicom
```

Get back to the Python interactive shell:

(Data source: <https://zenodo.org/record/16956#.XWc1b-P7Rpg>)

- a. Read the DICOM image `MR000008`, assign it to a variable `dcm`;
 - b. Use the print function to print `dcm`. What do you get?
 - c. Get the `pixel_array` from `dcm` and assign to `data`. Now `data` holds the imaging data.
 - d. Refer to the recording. Plot `data` using `matplotlib`;
 - e. Save `data` as `MR000008.tiff`;
 - f. Perform the same tasks in Question 1 on data and generate the corresponding files:
 - i. `MR000008_processed.tiff`;
 - ii. `MR000008_diff.tiff`;
 - iii. `MR000008_diff_rot.tiff`;
3. In the previous question, the `MR000008` is taken a 3D scan. The whole scan is in `SE000001.zip`. Unzip this file and inspect three more slices.